



Radar™

- ✓ **Non-Toxic**
- ✓ **Biodegradable**
- ✓ **Streak-Free**



INSTRUMENT & DEVICE ANTI-FOG PORTABLE SPRAY

Radar™ is a ready-to-use instrument and device anti-fog solution designed to ensure professionals see the clearest picture possible. Its nano-polymer formulation is designed to provide long-lasting anti-fogging on your screening devices.



FASTER



Convenient Sprayer

Pump sprayer is designed to make it quick and convenient to dispense the solution directly onto the surface.



Quick Fog Elimination

Immediately activates to create a barrier that resists fogging.



All Surfaces

Can be used on dental mirrors, loupes, visors, shields, clear masks, safety glasses, goggles and photography mirrors.

SAFER



Barrier

Forms a crystal clear barrier when used on intended surfaces.



Compatible

Can be used on glass, soft plastics, metals, silicone and acrylics and does not contain any chemicals which will break down plastics.



Environments

Effective in cold/dry or hot/humid environments.

KINDER



Environmentally Friendly

All products are made with biodegradable ingredients. Contains no toxic substances.



Recyclable

Packaging materials are fully recyclable.



Certified Biodegradable

All ingredients are USP Pharma grade and/or food grade quality.



INSTRUMENT & DEVICE

ANTI-FOG PORTABLE SPRAY

OUR MASCOT



Frogs are highly sensitive to their surroundings due to their permeable skin. The presence of a healthy frog population is a strong indicator of a balanced and healthy eco-system. Micrylium uses the frog as our mascot because we subscribe to creating products and packaging that are eco-friendly, biodegradable and sustainable.

SURFACE TENSION




Radar™ breaks down the surface tension of water, so instead of forming a fog of tiny droplets, the water spreads out into a thin, transparent sheet that doesn't obstruct vision. Radar™ also leaves behind a thin, invisible film or coating that repels water and prevents condensation from forming in the first place.

THE CONCEPT



Our goal is to provide solutions, not just to kill or inactivate, but to clean, protect and minimize the impact of regular sustained usage (people, equipment, professional and natural environments). Radar™ contains biodegradable surfactants, biological enzymes and naturally derived scents.

HOW TO USE



1 Spray over entire surface of mirror, loupes or other glass devices.



2 Wait 5-10 seconds and then wipe with LeCLOTH™ over surface.



3 This process functions in both warm and cold environments.



4 Proper application will result in streak-free visibility and a long lasting protection.



5 Use full strength, do not dilute.

PRODUCT SPECIFICATION DATA

| Item Number | Product Description | Packaging |
|-------------|-----------------------------------|------------|
| 05-FOGG-060 | Radar™ 60mL Mini Spray Bottle RTU | Case of 10 |

CUSTOMER REVIEWS


“I’ve been using Radar™ on my screening devices for the past few weeks and it’s made a huge difference. It dries quickly, doesn’t leave streaks, and gives me a consistently clear view. No more fogging during long procedures. Definitely a must-have for anyone who relies on precise visuals in their work.”

Karen W., Bright Smiles Family Dental, Ontario, Canada



RADCA3.31.10.25



| 1. IDENTIFICATION | | | | | | | | | |
|---|--|--|---------|--|--|---|--|---|--------------------------|
| Product Name | | Radar™ | | Manufacturer | | Micrylium Laboratories Inc. | | | |
| Registration | | CAN | Class I | Address | | 5000M Dufferin Street, Toronto, Canada, M3H 5T5 www.micrylium.com | | | |
| | | EU CE BAG | Class I | | | | | | |
| Indication | | Optical Medical Device Cleaner | | Phone | | 416-667-7040 | | | |
| Emergency Phone # | | CHEMTREC | | Fax | | 416-667-0071 | | | |
| | | | | 1-800-424-9300 | | CANUTEC | | 1-613-996-6666 | |
| 2. HAZARD IDENTIFICATION | | | | | | | | | |
| Symbol Pictogram | |  | | Signal Word | | Irritant | | | |
| | | | | Symbol | | Exclamation mark | | | |
| Classification | | Not Applicable | | | | | | | |
| Health Hazard | | No Serious Health Hazards | | | | Environmental Hazards | | Biodegradeable | |
| Precautionary & Hazard Statements | | <div><div><p>P102: Keep out of reach of children.</p><p>P301: IF SWALLOWED: Drink large quantities of water or milk.</p><p>P305: IF IN EYES: Flush eyes with large quantities of water.</p></div><div><p>H302: Harmful if swallowed.</p><p>H317: May cause an allergic skin reaction.</p><p>H320: Causes eye irritation.</p></div></div> | | | | | | | |
| 3. COMPOSITION | | | | | | | | | |
| Chemical | | | CAS # | | LD-50 (Oral, mg/kg) | | Concentration (%) | | |
| Ethanol | | | 64-17-5 | | 7,060 (Rat) | | 12.5% | | |
| 4. FIRST AID MEASURES | | | | | | | | | |
| Inhalation | | If breathing is difficult, remove individual to fresh air. | | | Ingestion | | Drink large quantities of milk or water to dilute. | | |
| Skin Contact | | No adverse effects. | | | Eye Contact | | Flush with plenty of water. | | |
| Most Important Symptoms and Effects (Acute and Delayed) | | | | | | | | | |
| May cause mild reversible irritation if there is eye contact or is inhaled. | | | | | | | | | |
| Indication of any Immediate Medical Attention and Special Treatment Needed | | | | | | | | | |
| Not Applicable. | | | | | | | | | |
| 5. FIREFIGHTING MEASURES | | | | | | | | | |
| Non-Flammable. | | | | | | | | | |
| 6. ACCIDENTAL RELEASE | | | | | | | | | |
| Use all means to prevent spillage. | | | | | | | | | |
| 7. HANDLING & STORAGE | | | | | | | | | |
| Store in a cool, dry, well-ventilated location. Keep away from heat, sparks and flames. DO NOT mix with bleach or peroxides. Storage and Transport: 5°- 30°C | | | | | | | | | |
| 8. EXPOSURE CONTROLS/ PERSONAL PROTECTION | | | | | | | | | |
| No specific measures required. | | | | | | | | | |
| 9. PHYSICAL AND CHEMICAL PROPERTIES | | | | | | | | | |
| Physical State | | Colour | Odour | Solidification point | Boiling point | Flash Point | Density g/ml @ 25°C | pH | Kinematic Viscosity@23°C |
| Transparent, Liquid | | Light Green | Apple | -5°C | 105°C | N/A | 0.99 | 6.5 | 12 mm²/s |
| 10. STABILITY AND REACTIVITY | | | | | | | | | |
| Stable under normal conditions. Incompatibility: Strong oxidants, acid chlorides, silver salts Decomposition: Products: CO ₂ , CO | | | | | | | | | |
| 11. TOXICOLOGICAL DATA | | | | | | | | | |
| Acute Dermal Toxicity | | LD ₅₀ >5000 mg/kg. Not found to be dermal sensitizer. | | | | Acute Oral | | N/A | |
| Ocular Irritation | | N/A | | | | Acute Inhalation Toxicity | | N/A | |
| Reproductive Hazards | | Ingestion/inhalation can be harmful. (TDLo 300mg/Kg Ethanol) | | | | Carcinogenicity | | Ingestion of Ethanol IARC Group 1. | |
| Tests Performed by Product Safety Labs, Dayton, NJ USA | | | | | | | | | |
| 12. ECOLOGICAL INFORMATION | | | | | | | | | |
| Surfactants are readily biodegradable in soil and water. Persistence unlikely based on available data. | | | | | | | | | |
| Ethanol | | EC50 (72h) = 275 mg/l (Chlorella vulgaris) | | Fathead minnow (Pimephales promelas) LC50 = 14200 Mg/L/96h | | Photobacterium Phosphoreum:EC50 = 34634 Mg/L/30 min Photobacterium Phosphoreum:EC50 = 35470 Mg/L/5 min | | EC50 = 9268 mg/L/48h EC50 = 10800 mg/L/24h | |
| 13. DISPOSAL CONSIDERATIONS | | | | | | | | | |
| Domestic. | | | | | | | | | |
| 14. TRANSPORT INFORMATION | | | | | | | | | |
| Not regulated | | | | | | | | | |
| N/A | | Land | | | Sea | | | Air (IATA) | |
| | | N/A | | | N/A | | | N/A | |
| | | N/A | | | N/A | | | N/A | |
| | | N/A | | | N/A | | | N/A | |
| 15. REGULATORY INFORMATION | | | | | | | | | |
| TSCA – No reporting required. | | | | | CERCLA – No hazardous pollutants or ozone depletion. | | | | |
| 16. OTHER INFORMATION | | | | | | | | | |
| The information and recommendations contained herein are based on information believed to be correct. It is offered in good faith, without guarantee. Micrylium Laboratories Inc. make no warranty expressed or implied. | | | | | | | | | |
| Effective Date: 2025/10/31 | | | | | Document: Radar 2.0 | | | | |